

Fully Automatic SUVR & Centiloid scale

R^2 with the Original Centiloid method : 0.98⁽¹⁾


Various A β tracer analysis

[¹¹C]PiB, [¹⁸F]Flutemetamol, [¹⁸F]Florbetaben, [¹⁸F]Florbetapir, [¹⁸F]NAV4694, [¹⁸F]Florapronol

3D Rendered Image


Customizable presets

FDG



- Specified for ROIs of Neurodegenerative disorders

General



- Accurate values for 100+ ROIs with or without MRI registration