Specification

Supported MRI data format

Image format	NIfTI - 1 (file extension : .nii)	
Туре	Structural MRI (T1-weighted)	
Slice thickness	Coronal ≤ 1.0 mm	
(Spacing)	Sagitta l ≤ 1.0 mm	
	Axia $I \le 1.6$ mm	
Fie l d Strength		
(Tes l a)	3.0 T	

Validated devices and protocols of MRI

Company	Model	Protocol
GE	Signa HDxt 1.5T	SPGR
Philips	Intera 1.5T	MPRAGE
	Intera 3T Ingenia 3T	TFE
Siemens	Skyra 3T Verio 3T	MPRAGE



NEUROPHET Inc. 12F, 124, Teheran-ro, Gangnam-gu, Seoul, Republic of Korea

T. +82 2 6954 7971 | F. +82 2 6954 7972 E. contact@neurophet.com

www.neurophet.com



Personalized tES Treatment tDCS/tACS Simulation Software tES LAB

neurophet

This version includes simulating electric field strength with stroke-affected region segmentation feature

Neurophet tES LAB is a complete, powerful tES simulation software


that allows you to simulate and analyze tES-induced electric fields in a personalized brain model using MRI.

Fully Automated Brain MRI Segmentation Including Stroke-affected Region



Deep learning-based Neurophet's technique segments skin, skull, CSF, white matter, grey matter, and even stroke-affected region using MRI.

Personalized Brain Modeling



An automated modeling algorithm generates a 3D brain model based on segmented tissue layers considering an individual's anatomical structures.

A Powerful tES Simulation Framework


The high-speed simulation feature calculates tDCS/tACS-induced electric field in the brain for advanced analysis and treatment planning.

- Easy, intuitive graphic user interface (GUI)
- The state-of-the-art visualization and analysis features (3D surface/cut-plane and ROI-based analysis)
- Electric field optimization features (single/multi-channel tES optimization)
- · Visual guidance for accurate electrode positioning

Additional Convenient Features

- Data export: segmentation label, 3D model, stimulation results
- Batch-processing for group data analysis

Fully Compatible with tES Parameters



- Customizable electrodes configuration (shape, size, number, and position of electrode)
- Electrode positioning system (10-20/10/5 EEG positioning system)

Neurophet also provides TMS LAB

for Personalized Transcranial Magnetic Stimulation Simulation and Planning

- Simulation & analysis for TMS-induced electric fields
- Optimization of a coil's position & orientation
- Interoperability with TMS Navigation*

*an upcoming function