


Personalized tES Treatment
tDCS/tACS Simulation
Software tES LAB

This version includes simulating electric field strength with stroke-affected region segmentation feature.

Neurophet tES LAB is a complete, powerful tES simulation software


that allows
tES-induced
personalized

Fully Automated Brain MRI Segmentation Including Stroke-affected Region




Deep learning-based Neurophet's technique segments skin, skull, CSF, white matter, grey matter, and even stroke-affected region using MRI.

Personalized Brain Modeling




An automated modeling algorithm generates a 3D brain model based on segmented tissue layers considering an individual's anatomical structures.

Fully Compatible with tES Parameters



- Customizable electrodes configuration (shape, size, number, and position of electrode)
- Electrode positioning system (10-20/10/5 EEG positioning system)



allows you to simulate and analyze
induced electric fields in a
personalized brain model using MRI.


A Powerful tES Simulation Framework

The high-speed simulation feature calculates
tDCS/tACS-induced electric field in the brain for
advanced analysis and treatment planning.

- Easy, intuitive graphic user interface (GUI)
- The state-of-the-art visualization and analysis features
(3D surface/cut-plane and ROI-based analysis)
- Electric field optimization features
(single/multi-channel tES optimization)
- Visual guidance for accurate electrode positioning

Additional Convenient Features

- Data export: segmentation label, 3D model, stimulation results
- Batch-processing for group data analysis



Neurophet also provides TMS LAB for Personalized Transcranial Magnetic Stimulation Simulation and Planning

- Simulation & analysis for TMS-induced electric fields
- Optimization of a coil's position & orientation
- Interoperability with TMS Navigation*

*an upcoming function

neurophet tES LAB

Specification

Supported MRI data format

Image format	NIFTI - 1 (file extension: .nii)
Type	Structural MRI (T1-weighted)
Slice thickness (Specimen)	Coronal ≤ 1.0 mm Sagittal ≤ 1.0 mm Axial ≤ 1.6 mm
Field Strength (Tesla)	1.5T 3.0T

Validated devices and protocols of MRI

Company	Model	Protocol
GE	Signa HDxt 1.5T	SPGR
Philips	Intera 1.5T	V/P-RAGE
	Intera 3T	TFE
	Ingenia 3T	
Siemens	Skyra 3T	V/P-RAGE
	Vero 3T	

neurophet

NEUROPHET Inc.

12F, 124, Teheran-ro, Gangnam-gu, Seoul, Republic of Korea

T. +82 2 6954 7971 | F. +82 2 6954 7972

E. contact@neurophet.com

www.neurophet.com

© 2023. NEUROPHET Inc., All Rights Reserved.

