



Fast & Accurate neurophet SCALE PET Analysis



High correlation
SUVR calculation in
multi PET tracer &
multimodal imaging
conditions⁴⁾



	MRI SEGMENTATION BY BRAIN REGIONS	REGISTRATION METHOD	FULL ANALYSIS TIME
CONVENTIONAL METHOD	8hours	Manual	8hrs or more
neurophet SCALE PET	1min	Fully Automatic	<div>PET+MR PET only</div> <div>PET+CT</div> <div>5min</div> <div>2min</div>

SCALE UP
Your Medical Insight



neurophet

NEUROPHET Inc.
12F, 124, Teheran-ro, Gangnam-gu, Seoul,
Republic of Korea

T. +82 2 6954 7971
F. +82 2 6954 7972
E. contact@neurophet.com

www.neurophet.com

©2024 NEUROPHET Inc., All rights reserved.

Visualization of PET Tracer Uptake
and SUVR Quantification

KOR


neurophet SCALE PET

MFDS(Korea) Registered

US-FDA 510(k) Cleared



PMDA(Japan) Registered

HSA(Singapore) Registered



neurophet

Amyloid



Fully Automatic SUVR & Centiloid scale

R² with the Original Centiloid method : 0.98⁴⁾


Various Aβ tracer analysis

[¹¹C]PiB, [¹⁸F]Flutemetamol, [¹⁸F]Florbetaben, [¹⁸F]Florbetapir, [¹⁸F]NAV4694, [¹⁸F]Florapronol

3D Rendered Image


Customizable presets

Tau




- SUVR for Tau PET staging (Braak stage ROIs & Temporal meta ROIs)
- [¹⁸F]Flortaucipir, [¹⁸F]MK-2640, [¹⁸F]PI-2620 etc

Dopamine




- Specified for striatum
- Asymmetry index
- Supportive for differential diagnosis of Parkinsonism(IPD, MSA, PSP) and ET

FDG



- Specified for ROIs of Neurodegenerative disorders


General



- Accurate values for 100+ ROIs with or without MRI registration


Neurophet SCALE PET supports clinical decision

Concordance of interpretation



Inter-rater Kappa on Amyloid PET positivity³⁾

*Validated in 3 university hospitals
*Compared with the standard reference determined with the agreement between two raters
***Rater 1:** Experienced Nuclear Medicine (≥10 years) / **Rater 2:** General Nuclear Medicine(≤3 years)



Ref. 1) Jiyeon Lee et al. Diagnostic 12(3), (2022) 623
2-3) Data on file, Neurophet Inc, CLIN-09-0001, CLIN-18-0001, GNT-04-0001, GNT-07-0001
4) Data on file, Neurophet Inc, 10620_DHF_UMKO